



Crucibles ♦ Tubes ♦ Plates ♦ More ...

Quartz Microscope Slides

Quartz microscope slides are used for high temperature applications up to 1200°C (2192°F) or microscopy applications which require UV transparency.

| Catalogue# | Length | | Width | | Thickness | | Price |
|------------|--------|------|-------|------|-----------|------|--------|
| | mm | inch | mm | inch | mm | inch | |
| FQ-S-001 | 25.4 | 1 | 25.4 | 1 | 1 | 0.04 | \$4.23 |
| FQ-S-002 | 50.8 | 2 | 25.4 | 1 | 1 | 0.04 | \$5.97 |
| FQ-S-003 | 76.2 | 3 | 25.4 | 1 | 1 | 0.04 | \$8.59 |

***For custom production & volume discount, please call (520) 514-1100
\$15.00 handling fees applied for orders less than \$100.00.
No minimum quantity requirements for standard products.***

AdValue Technology

Ceramic products for high-temperature and wear-resistant applications

Alumina Crucible • Alumina Boat • Alumina Tube • Alumina Plate

Quartz Crucible • Zirconia Crucible • Zirconia Boat • Custom Fabrication

[Http://www.advaluetech.com](http://www.advaluetech.com) • Tel: (520) 514-1100 • Fax: (520) 747-4024 • Email: sales@advaluetech.com